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Lipo3K® Transfection Reagent

Description

Lipo3K® is a newly developed and proprietary reagent for the transfection of nucleic acids into eukaryotic cells.

To obtain the highest transfection efficiency and low non-specific effects, optimize transfection conditions by varying DNA and Lipo3K® concentrations and cell number. Make DNA-Lipo3K complexes in serum-free medium such as Opti-MEM Reduced Serum Medium and add directly to cells in culture medium, in the presence or absence of serum/antibiotic. It is not necessary to remove complexes or change/add medium after transfection. The amount of Lipo3K Reagent for successful transfection varies. Start any new transfection by testing the recommended two concentrations of Lipo3K Reagent to determine an optimum amount.

Contents and Storage

Lipo3K® is supplied in liquid form in 2 bottles. Store at 4 ºC. DO NOT FREEZE.

Protocols : **Performing optimization first before transfection !**

To transfect cells with siRNA, follow the protocol as described for DNA.

1.Make sure that seed cells are greater than 70~90% confluent before transfection.

2.Dilute Lipo3K®-B Reagent in Opti-MEM Medium (2 tubes recommended) – Mix well

3. Prepare master mix of DNA by diluting DNA in Opti-MEM Medium, then add Lipo3K®-A Reagent – Mix well.

4. Add Diluted DNA to each tube of Diluted Lipo3K-B Reagent (1:1 ratio)

5. Incubate for 10–15 minutes at room temperature

6. Add DNA-lipid complex to cells

7. Incubate cells for 2–4 days at 37°C. Then, analyze transfected cells.

Procedure Details (Two Reaction Optimiz ation)			
Component	96-well	24-well	6-well
Adherent cells	1–4 × 10^4	0.5–2 × 10^5	0.25–1 × 10^6
Opti-MEM™ Medium	5 µL × 2	25 µL × 2	125 µL × 2
Lipo 3K-B Reagent	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL
Opti-MEM Medium	10 µL	50 µL	250 µL
DNA (0.5–5 µg/µL)	0.2 µg	1 µg	5 µg
Lipo 3K-A Reagent (2 µL/µg DNA)	0.4 µL	2 µL	10 µL
Diluted DNA (with Lipo3K-B Reagent)	5 µL	25 µL	125 µL
Diluted Lipo 3K-A Reagent	5 µL	25 µL	125 µL
Incubate for 10–15 minutes at room temperature			
component (per well)	96-well	24-well	6-well
DNA-lipid complex	10 µL	50 µL	250 µL
DNA amount	100 ng	500 ng	2500 ng
P3K™ Reagent	0.2 µL	1 µL	5 µL
Lipofectamine™ 3K Reagent used	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL
Incubate cells for 2–4 days at 37°C. Then, analyze transfected cells.			

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